

An Acoustical Approach to the Study of Marine Particle Dynamics

Clementina R. Russo

Dept. Physics & Astronomy

Emmanuel S. Boss

School of Marine Sciences

University of Maine

Many Thanks

Jim Loftin, Barbara Cole, Wayne H. Slade,
P.Hill, J. Newgard, J. Trowbridge
Office of Naval Research



Marine particle properties and dynamics

- Marine particles have optical and acoustical properties
 - Backscatter intensity is first-order correlated with mass concentration
 - Backscatter intensity per mass is dependent upon material properties of the scatterer
- Particle dynamics
 - Re-suspension/settling
 - Changes in particle size/composition
 - Aggregation
- Dynamic interactions may affect optical and acoustical properties of particles
 - Affected properties can affect predictability of suspended sediment concentration (implications for mass transport modeling)
 - Implications for applications in underwater communication



Quantifying marine particle properties and dynamics

- First-order variability in the water column
 - Previous studies have shown ABS, ADCP adept at measuring suspended sediment concentration
 - What about the ADV?
- Second-order variability in the water column
 - Optical properties of particles are affected by changes in composition and aggregation
 - What about the acoustic properties?
 - Disparate claims over instrument sensitivity to dynamical events


Laboratory Calibration of the ADV

- Assessed the backscattered acoustic intensity to increases in mass concentration
 - Concentration series performed with suspensions of particles of known size and composition
 - Multiple concentration experiments performed, each with a different size particle
- Characterized the mass-normalized backscatter to particle size
 - Compared results with acoustic scattering theory for suspensions of elastic particles with identical physical properties

Results, ADV Experiment



Results, ADV Experiment



ADV Stats

- Linear response to mass, up to 0.08gL^{-1}
- Mass-normalized backscatter relates to particle size as expected from theory, within uncertainties in K_t
- ADV shows promise in solving the inverse problem of determining concentrations

- **What about particle dynamics?**
 - Optics: aggregates alter PSDs, variable composition affects index of refraction
 - Disagreement about aggregation affecting acoustic properties

Aggregates in the Ocean

- Particulate composition

 - Organic

 - plant and animal matter

 - Inorganic


 - sand, mud, minerals, metals

- Particles aggregate

 - Following re-suspension/settling events

 - In expulsion from hydrothermal plumes


- Aggregates may be composed of various things, variously packaged



Aggregates in the sink

- Small particles of bentonite clay
 - Suspended in water, 0.04 g/L
 - Add salt, wait
 - Flocs form by differential sedimentation
- Measure properties over time
 - Acoustically – backscatter (ADV)
 - Optically – backscatter (BB3), size (LISST), optical attenuation (ac-9)
- Measure concentration over time
 - Suspended particle mass (SPM) pulled from sampling depth

Aggregation Experiment



BB3, optical backscatter




LISST particle sizer

ADV, acoustic backscatter




ac-9, optical attenuation meter

All sensors sampled from the same depth, SPM taken at depth, 9cm under the surface water

Results, Aggregation Experiment



Aggregation in the Field



Backscattering and Aggregation

- Results from lab experiment are consistent with field observations
- Aggregation affects optical and acoustic backscattering differently
- Mass-normalized acoustic backscattering is significantly decreased by aggregation
- Decrease may be attributed to acoustic absorption resulting from interactions between particles
 - Physical attributes of aggregates are different than solid particles
 - fractal, large fluid fractions promote different particle-particle/particle-fluid interactions
 - Implications for modeling acoustic backscattering
 - Implications for predicting suspended sediment concentrations

Acoustical Approach to Particle Dynamics

- Properties of particles are affected by dynamical interactions
 - Optical and acoustic backscattering are each affected differently
- Judicious co-deployments of optical and acoustical instrumentation can provide complimentary measurements
 - e.g., co-deployment may provide an indication of aggregation in-situ